



MultiAlign Tutorial 04 – Reviewing An Analysis

BRIAN LAMARCHE



About this tutorial

- This tutorial will describe the graphical user interface (GUI) after an analysis has been performed in great detail
- You will learn
 - Shortcuts for drilling through the data
 - How to interact with plots
 - Details about every view



Analysis overview

THIS SECTION PROVIDES A SET OF BACKGROUND INFORMATION TO HELP GUIDE YOU THROUGH THE REST OF THE TUTORIAL



Previous tutorials covered this material



Proudly Operated by Baffelle Since 1965

Wizard



1. Home Screen



2. Select Data



3. Set Parameters



4. Select Baseline and Mass Tag
Database



5. Set Analysis
Path and Name



Running Analysis



7. Analysis View Window






GUI Basics

UNDERSTANDING HOW TO INTERACT WITH PLOTS




Interacting with Plots



- This section describes how to interact with each plot.
- ► It is intended to provide basic to advanced functionality of the user interface
















Pacific Northwest





5500

















Context Menu (1)




- Right Click on the plot
- A context menu will appear
- Show All Data
 - Auto-zooms the plot to show all data points
- Copy All Data
 - Copies data from the plot to system clipboard for pasting into other software applications, e.g. Excel
- View Port
 - Allows you to step back to a previous zoom range

Pacific Northwest




Context Menu (2)



- Copy
 - Saves a WMF version of the plot, with the plot data still in raster form
- Copy Image
 - Will copy a rasterize image to the system clipboard
- Save Image
 - Saves the plot as an image to file
- Select All Series
 - Selects all series of points
- Selected Series
 - Details what series of points are selected



Measuring Distances



- Some plots will allow you to measure the distance between points
 - e.g. ppm and NET differences
- ▶ To Measure:
 - Hold Shift
 - Left Click
 - Drag
 - Release left mouse button
- ► The distance in the left window shows:
 - 1.86 ppm mass difference
 - .04 NET difference






Reviewing An Analysis


GLOBAL VIEW OF ANALYSIS

Open the analysis












Review of output




 Alignment related plots must exist in the Plots subfolder.

The analysis database file should also exist.

This is important later on because MultiAlign will search for these plots to create visualization elements in the GUI






- This is the main window that should appear
- If it does not, click on the analysis button at the bottom of the screen.





Analysis View – Major Views



If a mass tag database is used in the analysis, this tab allows investigation of each mass tag








Investigating clusters


THE CLUSTER TAB AND SELECTED CLUSTER VIEWS





- The clusters scatter plot window displays all of the clusters found in the analysis.
- Each blue dot represents a feature found across multiple datasets.
- This plot is interactive






- The clusters data grid shows all detailed information about each cluster.
- The columns are sort-able by clicking on the column header
- You can also reorder the column orders by clicking on a column header and dragging to its new location
- Left clicking on a row in the cluster data grid will display the details about the cluster below


Selected Cluster View Basics

Proudly Operated by Baffelle Since 1965







Tags and Clusters (1)













Tags and Clusters (3) – Shapes and Colors – Cluster Charts



Proudly Operated by Baffelle Since 1965

Shape Charge

+1

+2

+3

+4

X +5

Color Based on Type of Feature

LC-MS Features from selected cluster

LC-MS Features from nearby cluster (hollow)

Selected Cluster Centroid

Nearby Cluster Centroid

Features and Cluster Centroids



Matched mass tag




Unmatched mass tag nearby

Matched Mass Tag

LC-MS Features from nearby cluster (hollow diamond) – nearby mass tag unmatched

Mass Tags








Detailed information about each feature in cluster (consensus feature)

Number of MS
Features (i.e. number of parent scans)
feature appeared in


Number of MS/MS spectra associated with feature



LC-chromatogram of each charge state (depicted by color) for the selected feature



Features (2) - MS/MS






Abundance Profile

This slide currently shows the abundance (raw) from each dataset. Each dot represents the abundance for a feature from the dataset id below.


This plot is interactive

Future feature: organize each dataset by factor (e.g. age or time point)





MS/MS Spectra






Investigating mass tag matches

USING A MASS TAG CENTRIC APPROACH, FINDING HOW MANY CLUSTERS MATCH




Mass Tag View (1)



- The mass tag plot shows all mass tags present in the AMT Tag database used for alignment and/or peptide identification via STAC.
- Mass Tags are represented by orange diamonds.
- ► This plot is interactive




Mass Tag View (2)



- ► The mass tag histogram shows how many features (e.g. a peptide from an individual dataset) to the mass tag.
- ▶ This plot is interactive


Mass Tag View (3)



- ► The mass tag data grid shows detailed information about each mass tag.
- Left clicking on a row in the mass tag data grid displaysdetailed information about the
 - mass tag
















Pacific Northwest



Mass Tag Synopsis





Mass Tag - Tags and Clusters (1)



Proudly Operated by Battelle Since 1965


Clusters that match to this mass tag

And on other tab, nearby mass tags

Scatter plot that shows all matching clusters (red) and their features that comprise the cluster (blue)

Selected mass tag is displayed as orange diamonds. Nearby mass tags are displayed as hollow diamonds

NOTE: Feature and Cluster shapes are based on charge





Mass Tag - Tags and Clusters (2)



Proudly Operated by Baffelle Since 1965

Shape Charge

+1

+2

+3

+4

X +5

Color Based on Type of Feature

LC-MS Features of matched clusters

LC-MS Features from nearby cluster (hollow)

Matched Cluster Centroid

Nearby Unmatched Cluster Centroid

Features and Cluster Centroids




Selected mass tag



Nearby mass tags

Selected Mass Tag


LC-MS Features from nearby cluster (hollow diamond) – nearby mass tag unmatched Mass Tags





Mass Tag – Matched Proteins

Details about the proteins that this mass tag rolls into. Shown here is multiple mass tags because the mass tag (peptide in this case) selected has multiple conformations.





Proudly Operated by Baffelle Since 1965

Global Statistics Plots

UNDERSTANDING THE GLOBAL STATISTICS



Global Statistics Plots (no AMT tag database)



Proudly Operated by Baffelle Since 1965

Histogram of the number of datasets represented by each cluster

Histogram of the cluster sizes (number of features per cluster)



Histogram of charge states




Global Statistics Plots (no AMT tag database)



Proudly Operated by Battelle Since 1965

Histogram of the number of datasets represented by each cluster

Histogram of the cluster sizes (number of features per cluster)



Histogram of charge states

Histogram of mass tags that match to individual features (sum of cluster sizes for clusters that match to a mass tag)



Proudly Operated by Baffelle Since 1965

Datasets Plot

DATASET INFORMATION






Proudly Operated by Baffelle Since 1965

Allows the user to toggle between dataset plot view and data grid mode

MultiAlign







Dataset Plot View



Proudly Operated by Battelle Since 1965


These plots are created from the high resolution images stored in the Plots directory where the database was created.

Future feature will be to make these plots interactive





Dataset Plot View





Proudly Operated by **Battelle** Since 1965


Advanced Tab

APPLICATION CONFIGURATION



Advanced View





- This view allows you to customize the application and visualization defaults.
- This view currently has a minimal set of features displayed.



For more information see the MultiAlign website:

http://omics.pnl.gov/software/MultiAlign.php